

Refraction Stand
E-20 SE



**The modular design of the E20 SE allows you to design your own unit.
From the composition of the components to the colour and finish.**



Robust

A robust and solid metal frame is at the core of all bon refraction units. Made in Germany, bon units are built to endure the test of time.

- A sound investment in your future.



Flexible

The double sliding table will easily accommodate a larger instrument such as an autorefractor in the second position.

- Optimizing space and work flows where required.



Patient comfort.

The unique magnetic lock allows you to lock the table in any position allowing you to examine even the more corpulent patients without any embarrassing mishaps. Likewise the optional retractable arm rest can be moved up making the chair more accommodating.

- Putting you and your patient at ease.



Easy access for elderly patients.

The gently sloping base plate is designed to ensure easy access for patients. There is no step to trip over. The chair has been specifically designed so that your patient can sit down and get up ergonomically and easily without a struggle.

- Patient comfort guaranteed.



Your comfort.

The small details often make the difference. We have calculated exactly the turn on the swing over table to ensure you can achieve a smooth movement easily, without having to adjust your working position. All controls are logically placed and easily reached without you having to stretch or reach across your patient.

- Optimal working position.



Intelligent Design

Most refraction units require a safety off switch to prevent the chair from travelling too high and trapping the patient's legs. We put a little more thought into this aspect to ensure that this will never happen. No safety off switch is required. Likewise the chair travel has been so designed to ensure that you can test children and smaller adults without the need for an extra cushion.

- Patient comfort and safety is our common priority.



Logical Storage options

The E20 allows you to include an optional trial lens tray with integrated contact lens mirror. This is neatly stowed underneath the sliding table and can be swung out when required, allowing you a neat storage solution for trial lenses. You can practice lens insertion with your patient without the need for them to move.

- Saving time and space.

Computer friendly

In a digital world more and more computer equipment is entering the consulting room.

With the E20 you can make life easier by selecting the computer package. This is a neat solution comprised of;




A key board holder.

This is located under the sliding table and can be swung out when required. *) Your keyboard is at your fingertips as and when you need it.

Comfortable working position.

A computer rack for shuttle pcs provides neat storage of a compact desk top PC.



An Adjustable Monitor Mount.

This allows you to mount the monitor on the projector column. The flexible arm enables you to position the monitor so that you can work comfortably and share information and images with your patients.

*) This option cannot be fitted with the trial lens tray.

Design Solutions.



The projector column and phoropter arm can be height adapted to accommodate your phoropter of choice manual or automatic. Special adaptors which are colour matched and in keeping with the unit design enable us to fit any phoropter or projector.




The control panel is logically laid out and within easy reach. The intuitive lay out and the back ground illumination ensure that you will always find the correct switch even in a low light environment.

Chair up - Chair down - Spot light - Second instrument - Pulse for room lights - Auxiliary



The spot light is ideally positioned to provide ample light where your patient needs it, on the near vision card.

Compact



The E20 SE requires a maximum width of 1,50 m and is among the most compact of units available. And because the modular design allows you to incorporate all your extras, such as pc key board or trial set it really allows you to optimize your space. So if your room space is limited but you would like a versatile flexible unit that also looks good then you need look no further. - **E20 SE small wonder**



Customized Look

The control console is hand crafted from wood and with its curvaceous design creates a subtle look in an otherwise unsettling environment. We have tried really hard to create a feeling of a stylish piece of furniture to set your patients at ease.



Bring your own personal taste or corporate brand into the design. Select any colour combination from the RAL chart and match the upholstery from our extensive selection of samples.

- Personal branding.



Room planning service and installation.

With our wealth of experience we can assist you from the outset with your room planning and design. Our in house team can liaise directly with your architects to optimize space layout and work flows.

Whether you are relocating or redesigning our technical team will hold your hand and support you from A-Z. If we can assist we would be happy to - Let us know.



Footprint

