

## Manual Lensmeter **bon LM-2**

*This sturdy designed lensmeter is a well proven instrument for daily use. Its silver-grey design makes the bon LM-2 look very stylish.*



This device is suitable for spectacles measurement detection departments, spectacles processing factories, spectacles sales departments, departments of ophthalmology in hospitals and optical element factories to assay spherical lens diopter of ophthalmic lenses or lens, diopter, astigmatism axis of cylindric lens, diopter of prism and diopter of cornea-contact lens.

- measurement +/- 25D
- wide field of observation
- internal reading
- inclinable ocular
- including prism compensator



Specifications are subject to change without notice.