

## Pupillary Distance Meter

# PD-2



**Easy to operate.  
Precise and reliable results.**

- Uses energy saving, long lasting LED instead of common filament bulbs
- Range of Measurement:
  - Total Pupillary distance: 45 to 82 mm
  - Pupillary distance of left of right eye: 22.5 to 41 mm
- Measurement tolerance: less than 0.5 mm
- Automatic shut off: 1 minute after operation stops
- Voltage: 6V DC
- Size: 212 mm x 158 mm x 58 mm
- Weight: 0.6 kg
- Working distance measuring grades:
  - 30 cm, 35 cm, 40cm, 50 cm, 65 cm, 1 m, 2 m, infinite
- With fold out breathe shield
- Including batteries

Specifications, design and equipment are subject to change without any notice or obligation on the part of the manufacturer.