

Standard Instrument Tables



- ◁ **T-24**
Instrument table with electrically driven height adjustment 76 - 96 cm, central column, table top 48 x 62 cm



- ◁ **T-23**
Instrument table details as T-24 but with lateral column

High Quality - Selected materials and intelligent design are united to provide a secure and practical instrument table. The table base is designed with stability in mind. The castors are fitted with a break. The table top is edged with a black PVC coating for extra protection.

Powerful - The bon Instrument table with it's capacity to lift 200kg is truly versatile. It can be used for practically any instrument, even those heavier devices.
Voltage for all bon tables: 230V/50Hz

Quiet - The special housing and quiet motor result in a virtually noise free operation.

Elegant - The modern and simple design fits well into any environment.

Mobile - Equipped with concealed castors for easy movement.



- ◁ **T-24 V**
Instrument table with V shape table top for 2 large instruments.



- ◁ **TE-24**
Instrument table with double sliding table top for two large instruments table top 41 x 85 cm, height adjustment 81 - 101 cm