


bon T-50 instrument table



An extra wide lift range allows sitting or even standing working position.

Advantages at a glance:

- notably quiet, electro motorised height adjustment
- lift range from 71 to 120 cm
- suitable for exceptional heavy instruments weighing up to 150 kg
- two 240V power plugs below the table top
- colours: base: anthracite
tabletop: light grey (RAL 7035)
- symmetrical tabletop,
dimensions: 49 cm depth,
width: 51 cm (narrow range),
resp. 71 cm (broad distance)



Hint: optional larger table top available (ca. 90 x 49 cm)

Designs and details can be changed without prior notice.