

Ophthalmic Units

bon E-30 SE & bon E-40



bon has been making refraction furniture for two generations.

When you purchase a bon you are selecting over a quarter of a century of high quality engineering, a wealth of experience in practice lay out and a broad spectrum of innovative design ideas.

Ideas that can help you enhance your image and that of your practice.

Overview *bon E-30 SE & bon E-40*



Common Features as standard:

- Integrated desk unit.
- Double sliding table.
- Electromotorised height adjustable patient chair.
- Manual forwards backwards movement.
- One retractable armrest.
- Integral key pad control.
- Intelligent electronics
- Control of room lights via the unit.

◁ Both units provide a convenient service door that can be hinged down easily for comfortable service access.

Options common to both units:

- Special wood finishes.
- Phoropter arm.
- Projector column.
- A choice of cupboards to go beside or under the unit.
- Near reading spot light.
- ◁ Switchable holder for handheld instruments

Your Choice of Colours

Your style and your taste will be incorporated in your new unit. Select from over 250 RAL colours and numerous real wood veneers to match your practice profile. Upholstery of chairs and stools can be matched, too, choosing from 75 different colours. Besides that of course also special requests are possible.



Design

Design is a cornerstone of our philosophy. Our unique design provides you with a desk that is also a refraction unit. It is practical and it looks good. The open table design is practical for cleaning and for servicing.



Ergonomics

Ergonomics play a central role in any work environment. Your comfort and that of your patients are an essential given. The turn of a table, the location of controls, space for the patients spectacles or handbag; the ease of sitting down or standing up; the electrical elevation of the whole unit for correct working position.



Reliable

All our components are tested for high performance, because you need to be sure that your unit will perform time after time. With bon, reliability you can depend on is guaranteed.

Supported by our national service team of 8 highly qualified engineers you can be secure in the knowledge that be it installation or assistance, bon will be there to hold your hand.



Robust

At the core of all bon units is a rock solid metal mainframe, which ensures that they not only look good but they guarantee a top performance over and over again. Strong and solid, your bon unit is built to last.



Intelligent

Soft start electronics, programmable automatic lighting settings, optional programmable height adjustment, easy access for cleaning, infra red remote control, intelligent cabling are just some of the features that are designed in to make your life easier and your work day more enjoyable.

Ophthalmic Unit *bon E-30 SE*

Standard features:

- Electromotorised sliding table for two instruments.
- Electromotorised phoropter track.
- Infra red remote control.
- Beechwood finish


Options:

- Swing over trial lens tray.
- Optional electro-motorised height adjustment.
- Swing over disabled chair.
- A wide choice of wood & veneers and RAL finishes.



The E-30 desk curvacoues, stylish, practical. ▶

◀ Easy access for cleaning.



◀ Practical space for keys, mobile phone etc.

Easy to stand up and sit down. Especially helpful for elderly patients. ▶



Ophthalmic Unit *bon E-40*

Standard features:

- Manual sliding table for two instruments.
- Manual Phoropter track.
- Integral Trial lens tray with insert.

Options:

- Additional Writing desk.
- A wide choice of and or RAL finishes.



The integral spot light is ideally located for near vision testing. ▷

- ◁ No cables.
All cables and wiring are concealed –
no mess no clutter.



The patented angled phoropter track ensures that the phoropter comes in at the optimal angle. ▷



- ◁ Forwards backwards movement and the retractable armrest are ideal for the more corpulent patient.






Room Planning & Service

We can assist you with room design and electrical layouts. Our experienced team can speak directly with your architect. This seamless cooperation will optimize space utilization and patient flow management. Additionally our technical team is on hand to manage your installation and provide after sales service.

bon E-30SE Dimensions



bon E-40 Dimensions

