

## Pupillary Distance Meter PD-2

Easy to operate, precise and reliable results.




- Uses energy saving, long lasting LED instead of common filament bulbs
- Range of Measurement:  
Total Pupillary distance: 45 to 82 mm  
Pupillary distance of left of right eye: 22.5 to 41 mm
- Measurement tolerance: less than 0.5 mm
- Automatic shut off: 1 minute after operation stops
- Voltage: 6V DC
- Size: 212 mm x 158 mm x 58 mm
- Weight: 0.6 kg
- Working distance measuring grades:  
30 cm, 35 cm, 40cm, 50 cm, 65 cm, 1 m, 2 m, infinite
- With fold out breathe shield
- Including batteries

Specifications, design and equipment are subject to change without any notice or obligation on the part of the manufacturer.