


Refraction Stand

E-80

Your new workstation in unmatched quality



bon

The E-80 refraction unit

- a compact workstation in premium bon quality



Fascinating design

Design is an important matter for us when it comes to product development, especially when designing new, innovative examination units.

For instance, the E-80 is the "little sister" of the successful E-20SE design unit. This model is still the favorite workstation of hundreds of successful opticians and ophthalmologists.

The fresh appearance of the E-80 impresses with an attractive, pleasant color combination of light gray and sapphire blue. Its straightforward design adapts perfectly to any modern room and planning concept.

Fits all sizes

The seat height of the chair can be adjusted within a generous range.

The electric lifting column can handle a load of up to 140 kg. You can easily adjust the optimal examination height for the tabletop devices as well as for the phoropter. This allows you to examine large and small people equally comfortably - without an additional seat cushion. The seat cover is made of durable synthetic leather and is easy to clean.

Easy access

The E-80 deliberately has no bulky base plate. This clever design ensures barrier-free access without tripping hazards, even for older persons.



Flexible solutions

Any phoropter model can be adapted to the practical telescopic sliding arm, whether automatic or manual.

As an option to the telescopic sliding arm, the E-80 can also be equipped with a classic, swivelling AO-type counterbalanced arm.

The swivelling sliding table offers enough space for two large ophthalmic devices of your choice. Devices from other manufacturers can also be mounted without any problems.



Simply straightforward

The clearly arranged control panel of the E-80 is clearly limited to the really important functions: On-Off, Chair Up, Chair Down and an additional key for free assignment.

The illuminated buttons are easy to reach and extremely durable.



Reliable and durable


Like all bon units, the compact E-80 is manufactured solely in Germany to ensure durability and quality. A solid and extremely robust metal frame forms the core of each of our units.

The unit has integrated electronics that can automatically control the power supply to the table-top units, depending on which unit is currently in working position.

All instrument cables are invisibly installed inside the table. Sturdy cable guides safely protect them from bending and fraying.

Our certification as a medical device manufacturer in accordance with DIN EN ISO 13485 guarantees the fulfilment of the highest standards of quality and safety.


The bon E-80 fulfils the requirements of 2017/745 (EU) with flying colours.



A "bon" is always a wise investment in the future.

Space saving masterpiece

The bon E-80 requires a room width of just 1.43 m. That makes it one of the most compact units available and it fits comfortably even in small refraction rooms.



Room planning and installation

With our wealth of experience, we are happy to assist you with your planning from the very beginning. Our team can speak directly with your architects to work together to create the optimal space layout for effective workflow.

Whether you are relocating or redesigning your practice, our service team is happy to assist you from A to Z. Just ask us, we will be happy to support you!

Sit comfortably!

The "Samba 500" from the Danish manufacturer VELA Medical is the perfect companion for any refraction unit.

Its seat is exceptionally thickly padded.

So you sit extremely comfortably all day long.

The swivel stool has sturdy smooth-running castors and a practical height adjustment ring for easy access from all sides. The seat cover is made of easy-care and very robust imitation leather. The standard color is black. Many other colors are available on request.



E-80: bon Optic Vertriebsgesellschaft mbH, Stellmacherstr. 14, 23556 Lübeck
Samba 500: VELA Medical, Gøteborgvej 8-12, 9200 Aalborg SV, Denmark

Rev. 2022-09-02

bon Optic Vertriebsgesellschaft mbH

Stellmacherstr. 14 · 23556 Lübeck, Germany · Phone +49 451 80 9000
Fax +49 451 80 900-10 · E-Mail call@bon.de · Internet www.bon.de

The logo consists of the lowercase letters "bon" in a bold, blue, sans-serif font. The letters are slightly overlapping, with "b" on top, "o" in the middle, and "n" at the bottom.