

Instrument Tables
Series T-70



T-70
with symmetric table top, suitable for 1 large ophthalmic instrument, picture shows optional castors



T-70
with symmetric table top, suitable for 2 large ophthalmic instruments



T-70
with asymmetric table top, suitable for 1 large ophthalmic instrument



T-70
with special shaped table top for 1 large ophthalmic instrument plus **iMac** Computer, ideal for fundus camera **Cobra** or Scheimpflug camera **Sirius**

Instrument Tables Series T-70

The fresh colour design of the T-70 table series matches every interior perfectly.

table base: RAL 7024 graphite grey
table top: RAL 9003 signal white
cover plates: RAL 9007 grey aluminium

The pleasantly shaped table tops are available in different sizes, each suitable for your needs.

- Whether it's for one or two devices, with or without a computer.

The motor lifting column is very quiet and reliable. Depending on the requirements it is placed centrally or asymmetrically.

It's enormous lifting power corresponds to a weight of **over 150 kg**.

The modern concept of the T-70 gives you the option to upgrade up to two practical cable guides **through the table** top very easily. Therefore the table top has prefabricated perforations.

So all annoying power and data cables disappear elegantly next to the device.

Schuko sockets to **power your devices** are located directly below the table top.

Braked **castors** are available as optional accessory.

Specifications

height: 72 cm bis 115 cm
adjustment range: 43 cm
lifting power: 1500 N (ca. 153 kg)
power supply: 230 V / 50 Hz
connected load: 450W / max 2,0A

Specifications, design and equipment are subject to change without any notice or obligation on the part of the manufacturer.

